

Introducing *AG INDUSTRIES* **CPAP Tube Cleaning Wand**

Now you can easily clean the inside of your CPAP tube without any hassle or mess! This Tube Cleaning Wand quickly snaps together so you can clean your CPAP tube each morning after use. The disposable and durable cleaning pads attach to the end of your wand and slide along the length of any standard tube. When you're done, simply snap the wand apart and return to the convenient storage bag – great for travel!

FEATURES:

- Assembles to 40" to clean most tube lengths
- Durable and disposable cleaning pads
- Snaps apart for easy storage and travel
- Convenient storage bag


SOLD AS:

1 Week starter kit

Part # **AG-CLN-WA**

30 Day pad replacement

Part # **AG-CLN-WAREP**



Please contact customer service or your sales representative for more information on this product.